


Product Datasheet

DNAH6 antibody (orb1313106)

Description	Mouse monoclonal antibody to DNAH6
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human DNHL1 (NP_775916) produced in E.coli.
Target	DNHL1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	46.6 kDa
Uniprot ID	Q9C0G6
NCBI	NM_173645
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...