


Product Datasheet

Bone marrow stromal cell antigen 1 (BST1) antibody (orb1313107)

Description	Mouse monoclonal antibody to Bone marrow stromal cell antigen 1 (BST1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 29-318 of human BST1 (NP_004325) produced in E.coli.
Target	BST1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	35.5 kDa
Uniprot ID	Q10588
NCBI	NM_004334
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...