


Product Datasheet

KDM3A / JHDM2A (KDM3A) antibody (orb1313113)

Description	Mouse monoclonal antibody to KDM3A / JHDM2A (KDM3A)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 850-1321 of human KDM3A (NP_060903) produced in E.coli.
Target	KDM3A
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	147.2 kDa
Uniprot ID	Q9Y4C1
NCBI	NM_018433
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...