

Product Datasheet

Drebrin (DBN1) antibody (orb1313126)

Description	Mouse monoclonal antibody to Drebrin (DBN1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 2-240 of human DBN1 (NP_004386) produced in E.coli.
Target	DBN1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	71.3 kDa
Uniprot ID	Q16643
NCBI	NM_004395
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.

