


Product Datasheet

NCALD antibody (orb1313143)

Description	Mouse monoclonal antibody to NCALD
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human NCALD (NP_114430) produced in HEK293T cell.
Target	NCALD
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	22.1 kDa
Uniprot ID	P61601
NCBI	NM_032041
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...