


Product Datasheet

CREM antibody (orb1313147)

Description	Mouse monoclonal antibody to CREM
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human CREM (NP_874388) produced in E.coli.
Target	CREM
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	13.3 kDa
Uniprot ID	Q03060
NCBI	NM_182719
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...