


Product Datasheet

Apolipoprotein L 1 (APOL1) antibody (orb1313158)

Description	Mouse monoclonal antibody to Apolipoprotein L 1 (APOL1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 28-339 of human APOL1 (NP_003652) produced in E.coli.
Target	APOL1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	43.97 kDa
Uniprot ID	O14791
NCBI	NM_003661
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...