


Product Datasheet

ZNF8 antibody (orb1313182)

Description	Mouse monoclonal antibody to ZNF8
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-242 of human ZNF8 (NP_066575) produced in E.coli.
Target	ZNF8
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	64.8 kDa
Uniprot ID	P17098
NCBI	NM_021089
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...