


Product Datasheet

LAMP1 antibody (orb1313191)

Description	Mouse monoclonal antibody to LAMP1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 29-382 of human LAMP1 (NP_005552) produced in HEK293T cell.
Target	LAMP1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	44.88 kDa
Uniprot ID	P11279
NCBI	NM_005561
Dilution Range	WB 1:500, FLOW 1:25~100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...