


Product Datasheet

GIPC (GIPC1) antibody (orb1313197)

Description	Mouse monoclonal antibody to GIPC (GIPC1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human GIPC1 (NP_974199) produced in E.coli.
Target	GIPC1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	35.9 kDa
Uniprot ID	O14908
NCBI	NM_202470
Dilution Range	IHC 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...