


Product Datasheet

CRYGD antibody (orb1313200)

Description	Mouse monoclonal antibody to CRYGD
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human CRYGD (NP_008822) produced in HEK293T cell.
Target	CRYGD
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	20.6 kDa
Uniprot ID	P07320
NCBI	NM_006891
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...