


Product Datasheet

KAL1 (ANOS1) antibody (orb1313219)

Description	Mouse monoclonal antibody to KAL1 (ANOS1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 129-373 of human KAL1(NP_000207) produced in E.coli.
Target	ANOS1
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	76.11 kDa
Uniprot ID	P23352
NCBI	NM_000216
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...