


Product Datasheet

KIF6 antibody (orb1313275)

Description	Mouse monoclonal antibody to KIF6
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-350 of human KIF6 (NP_659464) produced in E.coli.
Target	KIF6
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	92.4 kDa
Uniprot ID	Q6ZMV9
NCBI	NM_145027
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...