

Product Datasheet

GPCR 2037 (GPR151) antibody (orb1313296)

Description Mouse monoclonal antibody to GPCR 2037 (GPR151)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen Full length human recombinant protein of human GPR151 (NP_919227) produced in HEK293T cell.

Target GPR151

Note For research use only

Isotype IgG1

Clonality Monoclonal


MW 46.5 kDa

Uniprot ID [Q8TDV0](#)

NCBI [NM_194251](#)

Dilution Range WB 1:500~2000

Expiration Date 12 months from date of receipt.



Western blot analysis of extracts (35ug)...