


Product Datasheet

ATP5ME antibody (orb1313297)

Description	Mouse monoclonal antibody to ATP5ME
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human ATP5I (NP_009031) produced in E.coli.
Target	ATP5ME
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	7.8 kDa
Uniprot ID	P56385
NCBI	NM_007100
Dilution Range	WB 1:500~2000, IHC 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...