


Product Datasheet

BPOZ (ABTB1) antibody (orb1313299)

Description	Mouse monoclonal antibody to BPOZ (ABTB1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human ABTB1 (NP_115937) produced in E.coli.
Target	ABTB1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	37.4 kDa
Uniprot ID	Q969K4
NCBI	NM_032548
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...