


Product Datasheet

MBD1 antibody (orb1313308)

Description	Mouse monoclonal antibody to MBD1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 222-556 of human MBD1 (NP_056723) produced in E.coli.
Target	MBD1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	61.1 kDa
Uniprot ID	Q9UIS9
NCBI	NM_015847
Dilution Range	WB 1:2000, IHC 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...