


Product Datasheet

Neurogranin (NRGN) antibody (orb1313319)

Description	Mouse monoclonal antibody to Neurogranin (NRGN)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-78 of human NRGN (NP_006167) produced in E.coli.
Target	NRGN
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	7.4 kDa
Uniprot ID	Q92686
NCBI	NM_006176
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...