



Product Datasheet

TSC22 domain family, member 4 (TSC22D4) antibody (orb1313333)

Description	Mouse monoclonal antibody to TSC22 domain family, member 4 (TSC22D4)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-395 of human TSC22D4 (NP_112197) produced in E.coli.
Target	TSC22D4
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	40.8 kDa
Uniprot ID	Q9Y3Q8
NCBI	NM_030935
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells u...



Western blot analysis of 2 cell lines us...