


Product Datasheet

PAX7 antibody (orb1313377)

Description	Mouse monoclonal antibody to PAX7
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 34-491 of human PAX7 (NP_039236) produced in E.coli.
Target	PAX7
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P23759
NCBI	NM_013945
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...