



Product Datasheet

TM4SF2 (TSPAN7) antibody (orb1313392)


Description	Mouse monoclonal antibody to TM4SF2 (TSPAN7)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human TSPAN7 (NP_004606) produced in HEK293T cell.
Target	TSPAN7
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P41732
NCBI	NM_004615
Dilution Range	WB 1:500, IHC 1:500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Carcinoma...



Immunohistochemical staining of Human li...



Immunohistochemical staining of Human en...