


Product Datasheet

MAGEB2 antibody (orb1313402)

Description	Mouse monoclonal antibody to MAGEB2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-224 of human MAGEB2 (NP_002355) produced in E.coli.
Target	MAGEB2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	35.1 kDa
Uniprot ID	O15479
NCBI	NM_002364
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...