


Product Datasheet

C2CD6 antibody (orb1313416)

Description	Mouse monoclonal antibody to C2CD6
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 307-623 of human ALS2CR11 (NP_689738) produced in E.coli.
Target	C2CD6
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	71 kDa
Uniprot ID	Q53TS8
NCBI	NM_152525
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...