


Product Datasheet

HOXA3 antibody (orb1313422)

Description	Mouse monoclonal antibody to HOXA3
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 40-276 of human HOXA3 (NP_109377) produced in E.coli.
Target	HOXA3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	46.2 kDa
Uniprot ID	O43365
NCBI	NM_030661
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...