


Product Datasheet

DTWD1 antibody (orb1313428)

Description	Mouse monoclonal antibody to DTWD1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human DTWD1 (NP_001138427) produced in E.coli.
Target	DTWD1
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	Q8N5C7
NCBI	NM_001144955
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...