


Product Datasheet

FOLR3 antibody (orb1313434)

Description	Mouse monoclonal antibody to FOLR3
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human FOLR3 (NP_000795) produced in E.coli.
Target	FOLR3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P41439
NCBI	NM_000804
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...