


Product Datasheet

TFEC antibody (orb1313487)

Description	Mouse monoclonal antibody to TFEC
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human TFEC (NP_001018068) produced in E.coli.
Target	TFEC
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	35.4 kDa
Uniprot ID	O14948
NCBI	NM_001018058
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...