


Product Datasheet

THAP8 antibody (orb1313511)

Description	Mouse monoclonal antibody to THAP8
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 86-274 of human THAP8 (NP_689871) produced in E.coli.
Target	THAP8
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q8NA92
NCBI	NM_152658
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...