


Product Datasheet

SPANXB1 antibody (orb1313543)

Description	Mouse monoclonal antibody to SPANXB1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human SPANXB1 (NP_115850) produced in E.coli.
Target	SPANXB1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q9NS25
NCBI	NM_032461
Dilution Range	WB 1:500~2000, IHC 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...