

Product Datasheet

PYGM antibody (orb1313585)

Description	Mouse monoclonal antibody to PYGM
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 698-842 of human PYGM (NP_005600) produced in E.coli.
Target	PYGM
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	96.9 kDa
Uniprot ID	P11217
NCBI	NM_005609
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...