


Product Datasheet

SorLA (SORL1) antibody (orb1313626)

Description	Mouse monoclonal antibody to SorLA (SORL1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 573-753 of human SORL1 (NP_003096) produced in E.coli.
Target	SORL1
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	Q92673
NCBI	NM_003105
Dilution Range	WB 1:500~2000, IHC 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...