


Product Datasheet

MSL2L1 (MSL2) antibody (orb1313631)

Description	Mouse monoclonal antibody to MSL2L1 (MSL2)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human MSL2 (NP_060603) produced in E.coli.
Target	MSL2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	62.4 kDa
Uniprot ID	Q9HCI7
NCBI	NM_018133
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...