


Product Datasheet

RAD21 antibody (orb1313633)

Description	Mouse monoclonal antibody to RAD21
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 193-458 of human RAD21 (NP_006256) produced in E.coli.
Target	RAD21
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	71.5 kDa
Uniprot ID	O60216
NCBI	NM_006265
Dilution Range	WB 1:2000, IHC 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...