


Product Datasheet

ABCF2 antibody (orb1313638)

Description	Mouse monoclonal antibody to ABCF2
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-235 of human ABCF2 (NP_009120) produced in E.coli.
Target	ABCF2
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	71.1 kDa
Uniprot ID	Q9UG63
NCBI	NM_007189
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...