


Product Datasheet

AKAP12 antibody (orb1313640)

Description	Mouse monoclonal antibody to AKAP12
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 266-557 of human AKAP12 (NP_005091) produced in E.coli.
Target	AKAP12
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	191.3 kDa
Uniprot ID	Q02952
NCBI	NM_005100
Dilution Range	WB 1:500~2000, IHC 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...