

Product Datasheet

TBX3 antibody (orb1313647)

Description	Mouse monoclonal antibody to TBX3
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 244-481 of human TBX3 (NP_057653) produced in E.coli.
Target	TBX3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	79.2 kDa
Uniprot ID	O15119
NCBI	NM_016569
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...