


Product Datasheet

CCDC27 antibody (orb1313658)

Description	Mouse monoclonal antibody to CCDC27
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human CCDC27 (NP_689705) produced in E.coli.
Target	CCDC27
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	75.2 kDa
Uniprot ID	Q2M243
NCBI	NM_152492
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
HEK293T
cells w...