

AGO1 antibody (orb1313671)

Western blot analysis showing protein levels of Hsp70, Hsp90, and p70S6 in Hep2, Hela, M7, A549, C2C12, A431, MDCK, PC12, and MCF cell lines. Molecular weight markers (191, 97, 66, 51, 39, 28, 19, 14 kDa) are indicated on the left. Hsp70 and Hsp90 bands are visible in all lanes, while p70S6 bands are visible in all lanes except MCF.

Western blot
analysis of
HEK293T
cells w...