


Product Datasheet

CIB1 antibody (orb1313690)

Description	Mouse monoclonal antibody to CIB1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 2-191 of human CIB1 (NP_006375) produced in E.coli.
Target	CIB1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	21.5 kDa
Uniprot ID	Q99828
NCBI	NM_006384
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...