


Product Datasheet

IL17B antibody (orb1313702)

Description	Mouse monoclonal antibody to IL17B
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human IL17B (NP_055258) produced in E.coli.
Target	IL17B
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	20.44 kDa
Uniprot ID	Q9UHF5
NCBI	NM_014443
Dilution Range	WB 1:2000, IHC 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...