


Product Datasheet

ZNF670 antibody (orb1313704)

Description	Mouse monoclonal antibody to ZNF670
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human ZNF670 (NP_149990) produced in E.coli.
Target	ZNF670
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	44.4 kDa
Uniprot ID	Q9BS34
NCBI	NM_033213
Dilution Range	WB 1:500~2000, IHC 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...