

Product Datasheet

CSPS (SULT1A3) antibody (orb1313707)

Description Mouse monoclonal antibody to CSPS (SULT1A3)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications IHC

Immunogen Full length human recombinant protein of human SULT1A3 (NP_003157) produced in E.coli.

Target SULT1A3

Note For research use only

Isotype IgG1

Clonality Monoclonal


MW 34 kDa

Uniprot ID [P50224](#)

NCBI [NM_003166](#)

Dilution Range IHC 1:500

Expiration Date 12 months from date of receipt.



Immunohistochemical staining of paraffin...