

Product Datasheet

CSPS (SULT1A3) antibody (orb1313708)

Description Mouse monoclonal antibody to CSPS (SULT1A3)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications IHC, WB

Immunogen Full length human recombinant protein of human SULT1A3 (NP_003157) produced in E.coli.

Target SULT1A3

Note For research use only

Isotype IgG1

Clonality Monoclonal


MW 34 kDa

Uniprot ID [P50224](#)

NCBI [NM_003166](#)

Dilution Range WB 1:500~2000, IHC 1:500

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...