


Product Datasheet

MMGT1 antibody (orb1313710)

Description	Mouse monoclonal antibody to MMGT1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human MMGT1 (NP_775741) produced in E.coli.
Target	MMGT1
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	14.5 kDa
Uniprot ID	Q8N4V1
NCBI	NM_173470
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...