


Product Datasheet

CEP72 antibody (orb1313717)

Description	Mouse monoclonal antibody to CEP72
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human CEP72 (NP_060610) produced in E.coli.
Target	CEP72
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	71.5 kDa
Uniprot ID	Q9P209
NCBI	NM_018140
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...