


Product Datasheet

DUSP5 antibody (orb1313780)

Description	Mouse monoclonal antibody to DUSP5
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 125-384 of human DUSP5 (NP_004410) produced in E.coli.
Target	DUSP5
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	41.9 kDa
Uniprot ID	Q16690
NCBI	NM_004419
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...