


Product Datasheet

ASH2L antibody (orb1313796)

Description	Mouse monoclonal antibody to ASH2L
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 347-628 of human ASH2L (NP_004665) produced in E.coli.
Target	ASH2L
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	68.5 kDa
Uniprot ID	Q9UBL3
NCBI	NM_004674
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...