


Product Datasheet

IMPDH2 antibody (orb1313806)

Description	Mouse monoclonal antibody to IMPDH2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human IMPDH2 (NP_000875) produced in HEK293T cell.
Target	IMPDH2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	55.6 kDa
Uniprot ID	P12268
NCBI	NM_000884
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...