


Product Datasheet

PTF1A antibody (orb1313826)

Description	Mouse monoclonal antibody to PTF1A
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human PTF1A (NP_835455) produced in E.coli.
Target	PTF1A
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	34.8 kDa
Uniprot ID	Q7RTS3
NCBI	NM_178161
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...