


Product Datasheet

SLU7 antibody (orb1313836)

Description	Mouse monoclonal antibody to SLU7
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-300 of human SLU7 (NP_006416) produced in E.coli.
Target	SLU7
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	68.2 kDa
Uniprot ID	O95391
NCBI	NM_006425
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...